


# Dornhoffer



Grace Medical  
ENT. Redefined.

## Implant Line



HA contoured head design

Designed to lock between the stapes and malleus

Available with windows for increased visualization

Malleus groove

Footplate shoe, adjustable to fit between stapes crus remnants

636-080 | .8mm Footplate Shoe

Ref. #	Type	Ref. #	Type
635	Adjustable Total	637	Adjustable Total w/ Windows
695	Adjustable Partial	687	Adjustable Partial w/ Windows
690	2mm Fixed Length Partial	691	2mm Fixed Length Partial w/ Windows

626 | Adjustable Total w/ Windows and Integrated Footplate Shoe

Designed in conjunction with John Dornhoffer, MD, Little Rock, AR

8500 Wolf Lake Drive, Suite 110 • Memphis, TN 38133 • USA Toll Free 1.866.472.2363 • [www.gracemedical.com](http://www.gracemedical.com)

LIT0064 CID2892