

MELTI Piston




Large Loop For Easy Placement On Incus

2x Wider Than Standard Platinum Ribbon Prevents Rolling When Crimping

Smooth Transitioned Depth Marker Prevents Interference With Fenestration in Footplate

0.5mm Depth Marker To Aid Surgeon In Visualizing Placement Into the Vestibule



Ref. #	Materials	Shaft Dia.	Length
418-375	Platinum & Titanium	0.5mm	3.75mm
418-400	Platinum & Titanium	0.5mm	4.00mm
418-425	Platinum & Titanium	0.5mm	4.25mm
418-450	Platinum & Titanium	0.5mm	4.50mm
418-475	Platinum & Titanium	0.5mm	4.75mm

Designed in conjunction with Dennis Bojrab, MD and Michael LaRouere, MD, Detroit, MI

8500 Wolf Lake Drive, Suite 110 • Memphis, TN 38133 • USA Toll Free 1.866.472.2363 • www.gracemedical.com

LIT0096 CID2892